



This **Universal Rail Clamp** can position surgical devices anywhere along standard operating room table side rails. Use the Schure Socket to halt all movement on mounting posts as well.

## Features

- Clark Sockets by SchureMed are designed for a wide-range of purposes; affixing most operating room table accessories to most OR table designs and brands
- Spring-loaded and attaches to O.R. table side rails anywhere along the rail
- It accepts both 5/8" (1.6 cm) diameter and 1" x 3/8" (2.5 cm x .95 cm) flat bar mounting posts
- Holds all accessories in the vertical and horizontal position
- Patient weight capacity is 350 lbs (158 kg)
- Locks securely with a simple twist of the ergonomic handle

