



The **Valleylab CUSA 200** is an ultrasonic surgical system. The CUSA 200 is offered by Soma Tech Intl refurbished up to 50% below OEM prices with the same service and warranty contracts available. This ultrasonic system provides great tissue fragmentation with consistent performance.

## Features

- Provides quick and selective tissue fragmentation.
- Tactile feedback and exceptional visualization of the operative site.
- Great tip integrity that provides consistence performance and minimal potential for tissue damage.
- Integrated suction and irrigation that allows for flushing and aspiration without introducing other instruments into the operative site.



## Specifications

### Dimensions

**Height:** 76" (193 cm)  
**Width:** 26" (66 cm)  
**Depth:** 23" (58 cm)  
**Weight:** 226 lbs (102 kg)

### Ultrasonic System

**Frequency:** 23 kHz  
**Maximum Tip Amplitude:** Up to 355 microns  
**Oscillation System:** Magnetostrictive  
**CAVI-Pulse:** Four distinct power cycles

### Fluidic System

**Aspiration:** Up to 24 in. Hg vacuum  
**Irrigation Rate:** 1.5 to 50 cc/min  
Coaxial Irrigation.  
Automatic tissue release from aspiration

### Electrical

3 prong hospital grade plug 100/110/120 VAC  
15 amps  
50-60 Hz  
  
200/220/240 VAC  
10 amps  
50-60 Hz

### Handpieces and Tips

**PMT Straight Handpiece:** Length: 9" (23 cm)  
**PMT Angled Handpiece:** Length: 9" (23 cm)  
**PMT Extenders**  
**Straight:** Length: 4" (10.2 cm)  
**Curved:** Length: 4" (10.2 cm)  
**Titanium Tips**  
**Standard:** Length: 3.1" (7.9 cm) Outside Diameter: 0.1" (0.254 cm)  
**Micro Tip:** Length 3.6" (9.2 cm) Outside Diameter: 0.076" (0.193 cm)  
**Laparoscopic:** Length 12.3" (31.2 cm) Outside Diameter: 0.1" (0.254 cm)

