



The **Verathon GlideScope AVL** video laryngoscope system delivers the digital video technology and functionality physicians expect from advanced video laryngoscopy. This Laryngoscope is ideal for difficult and unpredictable airways, and provides a consistently clear view that enables quick intubation in unpredictable conditions.

Features

- Digital Color Monitor.
- Integrated real-time recording.
- Single-use blade.
- Ergonomic Design.
- Anti-fog technology.



Specifications

Dimensions

Height: 6.8 in (17.4 cm)

Width: 8.8 in (22.3 cm)

Depth: 3.1 in (8 cm)

Weight: 2.2 lbs (1 kg)

Display

Type: TFT Color VGA

Resolution: 640 x 480 px

Size: 6.4 in (16.3 cm)

Battery

Type: Lithium-ion

Battery Life: ~ 90 minutes

Charging Time: less than 6 hours

Capacity: 2150 mAh

Voltage: 7.2 V

Weight: 3.17 oz (90 g)