



GlideScope Go is a durable, portable, high-resolution video laryngoscope system featuring a scratch-resistant 3.5-inch color display with anti-glare coating and a wide field of view. Designed for improved visualization and rapid, first-pass intubation, it is compatible with GlideScope Spectrum single-use blades, for patients ranging from neonate to large adults.

### Convenience

Get the confidence of GlideScope wherever and whenever you need it, with the new portable GlideScope Go™. Compact, durable, easy to charge and ready on demand, GlideScope Go simplifies maintenance with a fully submersible screen and battery, and fully disposable, single-use blades. An essential complement to your GlideScope cart-based system, this handheld system is also ideal for use in small spaces, emergent procedures and whenever your environment demands mobility.

### Compatibility

Access the portfolio\* of GlideScope Spectrum fully disposable, single-use blades—each compatible with your GlideScope cart-based system—and generate clear airway views.

### Clarity

Get more of the light you want and less of light you don't. GlideScope Go delivers a high-resolution image powered by the latest advancements in lighting and camera technology—including Dynamic Light Control™ and Ambient Light Reduction™. With a glare-resistant, 3.5" screen and wide field of view, this intuitive system is designed for improved visualization and rapid first-pass success in patients ranging from neonates to large adults.





## Features

- Access the GlideScope Spectrum blade portfolio
- Fully disposable, single-use blades to support infection control initiatives
- 3.5" landscape, color display with vertical tilt adjustment for maximum flexibility
- Easy to use and easy to train, with on/off button and intuitive display
- Hot-swappable functionality enables blade changes without powering down
- Integrated battery charges to 80% in 2 hours and delivers a minimum of 100 minutes of continuous use on a full charge
- Automatic recording with removable micro-USB drive
- Scratch-resistant screen with anti-glare coating for reliable clarity in any environment
- Fully submersible video monitor and battery to simplify cleaning
- Dynamic Light Control™ and Ambient Light Reduction™ for enhanced airway views
- Operates on corded, external power should battery be discharged
- Remaining Charge indicator eliminates guessing
- Rugged design ensures functionality even when dropped from 4 feet
- IP67 rating ensures full dust protection and reliable operation in wet conditions
- Backed by a full warranty and customer support

