

Features

-
- A GlideScope video laryngoscope system. The main unit is blue and white, featuring a large screen and a control panel with buttons for MENU, UP, DOWN, ON/OFF, BATTERY INDICATOR, and CHARGE STATUS. A blue handle with a camera at the end is connected to the unit.

Specifications

Display Monitor

LCD TFT Panel: 320 x 240

Monitor: 6.4 in

Height: 167.3 mm

Width: 210.0 mm

Depth: 87.6 mm

Weight : 1.4 kg

GlideRite® Rigid Stylet

Handle Length: 82 mm

Handle Width: 16 mm

Stylet Rod Length: 266 mm

Distal Tip Diameter: 5 mm

Power

AC power cord: 4.5 m (15 ft)

Video Output Cable: (for connection to NTSC-compatible devices)
0.8 m (2.7 ft)

Video Cable: 1.8 m (6 ft)

