



The **Welch Allyn CP 150** is an EKG system that can also measure spirometry. The CP 150 comes with a 7-inch color touchscreen display with a keyboard. The Welch Allyn system can easily provide full-page ECG readings with the touch of a button. With options for 3, 6, or 12-lead EKG options. The Resting electrocardiograph system can store up to 100 readings on its internal memory.

## Features

- 3, 6, and 12-lead EKG readings
- 7-inch color touchscreen display
- Light weight battery operated system
- Can store up to 100 readings
- Offers Spirometry measurements
- USB and Ethernet connectivity to manage patient data

