



Designed specifically for demanding viewing conditions, with the ability to control the light output, the Green Series 600 Minor Procedure Light features three white LEDs in a wider head with integrated handle for easy maneuverability. LEDs do not require replacement and are a brighter light than halogen, while consuming less energy for improved facility efficiencies. Now you can see true tissue color rendition while reducing your current cost of ownership. The Green Series Lights were designed to improve patient exams while reducing environmental impact.

## Features

- Three White LEDs for greater intensity than halogen
- Minimal heat dissipation
- Wider head, designed specifically for demanding viewing conditions
- Uniform spot and superior brightness for minor procedures and exams in any healthcare environment
- Touchless power controls and adjustable light output
- Disposable sheath available to help reduce the risk of infection due to cross-contamination
- Rugged, ergonomic design allows for improved performance and long-lasting durability
- Low cost of ownership due to elimination of bulb replacement and lower energy consumption
- Easy to assemble and clean
- 5 year Welch Allyn warranty to include the LED light source and power supply only



# Specifications

## Physical Dimensions

### Weight:

Total Table/Wall Mount: 6.45 lbs (2.93 kg)

Total Mobile Stand: 17.35 lbs (7.87 kg)

### Control Box:

Height: 10.83 in (27.5 cm)

Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm)

### Luminaire:

Height: 13.58 in (34.5 cm)

Width: 7.56 in (19.2 cm) Depth: 3.39 in (8.6 cm)

### Mobile Stand:

Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm)

**Total Wall Mount:** 33.27 in (84.5 cm)

**Total Mobile Stand:** 58.66 in (149.0 cm)

**Table/Wall and Mobile Stand to Luminaire Center:** 31.1 in (79 cm)

## Number of LEDs

3 White LEDs

## Color Temperature

5,500° K

## Light Intensity @ Typical Working Distance

50,000 lux (4,645 footcandles)

## LED Life

50,000 hours

## Working Distance

12 in (30.5 cm) to 20 in (50 cm)

## Spot Size at Working Distance

6 in (15 cm) to 7 in (18 cm)

## Electrical

Input: 100-240 V~50 Hz-60 Hz,  
0.65 A max

## Mounting Options

Table/Wall Mount  
Mobile Stand

## Control Box

Touchless On/Off  
Control with Dimming