

Welch Allyn Quinton TM55 Stress Treadmill



The **Welch Allyn Quinton TM55** Stress Treadmill where simplicity and precision come together seamlessly. With its whisper-quiet operation, tests are conducted in a serene environment, enhancing patient focus. Effortless operation is at your fingertips, eliminating complexities. At its core, the Ultradrive system ensures impeccable control, delivering accuracy even at slow speeds and high weights. Elevate your stress testing with the TM55 for unmatched ease and reliability.

Features

- Easy-to-read cardiac stress systems simplify the process of starting, controlling, and stopping the treadmill.
- Clear displays allow you to easily locate test duration, speed, and grade.
- The stress system allows you to bring your patients to a comfortable stop, gradually reducing speed and lowering the grade to zero.
- Optional rapid deceleration capability is available for the times you need to slow and stop the treadmill quickly, such as stress echo peak exercise imaging.



Welch Allyn Quinton TM55 Stress Treadmill

Specifications

Dimensions

Height: 46.5 in (118 cm)

Width: 29.9 in (75.9 cm)

Length: 80.2 in (204 cm)

Weight: 352 lb (160 kg)

User Weight Capacity

500 lb (227 kg)

Speed Range

0.8 mph to 9.6 mph (1.3 km/h to 15.4 km/h) in 0.1 mph increments

Incline Range

0% to 25%, in 0.5% increments

Power

Low Voltage — 100/120 VAC @ 20 A; 50/60 Hz; NEMA 5-20R non-locking receptacle; Single-phase dedicated circuit

High Voltage — 200/230 VAC @ 10 A; 50/60 Hz; NEMA 6-15R non-locking receptacle; Single-phase dedicated circuit

