



The **Zeiss OPMI 11 Surgical Microscope** is mounted on an S21 stand. The OPMI 11 is an ENT Surgical microscope and is capable to accommodate for simple tubes as well as more advanced procedures. These procedures include Ear, Nose and Throat. The Microscope offers a single, fiberoptic light source with a mercury switch for on/off. With a large 65 mm field -of-view and spot illumination, and a five step magnification (4x, 6x, 10x, 16x, and 25 x), fine focus and a larger diameter objective lens and wide field eyepiece.

Features

- Smooth Fine-Focus Adjustment
- Fiberoptic Illumination with Mercury Switch
- Durable Metal Components
- Compatible with all various mounts



Specifications

General

Binocular: Straight or Inclined; Inclinable Upgrade Available

Objective Lens: Your Choice: 150mm-400mm

Eyepieces: 10X Widefield; 12.5X Available

Pupillary Distance: 50mm - 75mm

Filter: Standard UV; Green Filter Available

Lightsource: 12V 100W

Illuminated Field: Up to 65mm

Power Source: 110-240V ; 50 - 60Hz

