



The **Zeiss OPMI CS NC-2 Surgical Microscope** is ideal for spinal surgeries. The CS NC-2 comes with two independent zooms and internal focusing systems for Face to Face optics. There's also an assistant scope that comes with a 10X wide field eyepieces and a Xenon light Source.

Features

- Maximum Reach and Maneuverability
- Opmi CS Microscope
- NC2 Contraves floorstand Optional Automatic Leveling System
- 300 w Xenon light source w/Halogen backup
- Fiber Optic Cable
- Variable Working Distance
- Face to Face Optics
- Variable Focus: 200mm – 400mm
- Handlebar Controls for Zoom, Focus and Brake Release
- 2 0-180 degree Inclinable binoculars with 10x wide field Eyepieces
- 1 Side Assistant's Scope with 10x widefield Eyepieces
- Zeiss Verioskop



Specifications

General

Binocular Inclinable: 0 - 180°; face to face packages available
Objective Lens: Integrated 207 mm - 400 mm working distance
Magnification: High resolution, motorized via handlebars
Focus: Micro fine, motorized via handlebars
Eyepieces: 10x widefield Pupillary Distance: 55 mm - 75 mm
Filter: Standard UV
Lightsource: 300w Xenon with halogen back-up
Illuminated Field: Constant at all magnification levels
Power Source: 110-240V ; 50 - 60Hz
Mounting System: Contraves
Horizontal Reach: Arm = 62" | End of optics = 69"
Vertical Range (Macro focus): 75"
Arm: Free-Floating via handlebars
Shipping Weight: ±600lb