



The **Zeiss OPMI LUMERA S7** features more clarity, contrast, and stability than previous models. The LUMERA i allows you see intricacies of the ocular anatomy that you have never seen before with the help from Stereo Coaxial Illumination (SCI). Contours appear razor sharp. This microscope allows you to easily recognize different structures and fluidics, even with dense cataracts. The LUMERA S7 has the ability to equip a stereo co-observation tube allowing a second person to see the surgical field at the same magnification level.

Features

- Targeted surgical area shines brightly and vividly.
- DeepView allows you to switch between maximum transmission and the maximum depth of field, whichever is required.
- S7 manual floor stand that is compact and maneuverable.
- Motorized zoom.
- Zoom ratio: 1:6.
- Apochromatic optics.
- Focusing Range: 50 mm.
- SCI: Red reflex illumination.



