



The Zeiss OPMI MDU is an ophthalmic surgical microscope. It can be used for Cataracts as well as other eye surgeries and procedures for hand and plastic surgery. It has a removable optical components that allow for a face to face bridge. The MDU delivers a very consistent red reflex and great Zeiss optical clarity. The S5 microscope is a compact, manually balanced stand.

Features

- Great Optical Clarity
- Usable for Multiple Disciplines
- Zeiss Opmi MDU
- Zeiss XY2
- Zeiss F170 0-180 degree inclinable binos
- Zeiss F170 Inclined binos
- 4 10x widefield eyepieces
- Zeiss beamsplitter
- Multi-axis stereo observer tube
- Zeiss objective lens 175mm
- Zeiss S5 floorstand
- Zeiss footswitch
- Sterilizable covers

