



The **Zeiss OPMI 6 Surgical Microscope** is mainly used for eye surgeries. The OPMI 6 Surgical Microscope has a sharp optical view, a terrific Depth of Focus, and a compact microscope body. It features inclined binoculars that come with eyepieces. The OPMI 6 has superb clarity and enhanced red-reflex. The Zeiss Microscope has the newer Visual series technology, to improve the microscopes simplicity and it's performance at the same time. This system is one the least expensive microscopes for eye surgery.

Features

- Great Value on a Full Ophthalmic Scope
- Reliable Red-Reflex
- Zeiss Opmi 6
- Zeiss XY2
- 0-180 degree inclinable binos
- Zeiss F125 inclined binos
- 4 12.5x widefield eyepieces
- Zeiss beamsplitter
- Stereo observer tube
- Zeiss objective lens 175mm
- Zeiss S3 floorstand
- Zeiss footswitch
- Sterilizable covers



Specifications

Microscope

Binocular: Inclined ; f=125mm | Other options: 0-60 Inclinable ; f=160mm or 0-60 Widefield ; f=170mm

Objective Lens: 175mm, 200mm

Magnification: Motorized 1:4, .5...2.0

Focus: Motorized with 50mm travel

Eyepieces: 12.5X HEP Push-In

Pupillary Distance: 50mm - 70mm

Filter: Standard UV

Illumination

Lightsource: External 15V | 150W

Illuminated Field: 48mm

Power Source: 110 - 240V ; 50-60Hz

Mounting System

Type: Standard 1880 Floorstand

Horizontal Reach: Arm = 24" ; End of Optics = 27"

Vertical Range: 72"

Rotation Angle: 350

Shipping Weight: 220 lbs with weights for counterbalancing