



**Zeiss 6S-FC microscope** has always been reliable both in optical clarity and with enhancing the red-reflex. The technology was applied to the newer Visu series because of its simplicity and performance. The Zeiss 6SFC and the S3 stand are ideal for an ophthalmic microscope.

## Features

- Great Value on a Full Ophthalmic Scope
- Reliable Red-Reflex
- Zeiss Opmi 6s
- Zeiss XY2
- 0-180 degree inclinable binos
- Zeiss F125 inclined binos
- 4 12.5x widefield eyepieces
- Zeiss beamsplitter
- Stereo observer tube
- Zeiss objective lens 175mm
- Zeiss S3 floorstand
- Zeiss footswitch
- Sterilizable covers

