

# Zimmer IntelliCart System Waste Management



The **Zimmer IntelliCart System** sets the industry standard with its impressive 34-liter fluid capacity, ultra-quiet vacuum pump, obstruction-free suction manifolds, and portable smoke evacuation capability. Experience the ease and convenience of effectively managing surgical fluid and smoke waste hazard.

## Features

- Redesigned Vacuum Pump that is extra quiet with HEPA air filtration, increased flow capabilities, and odor control.
- Dual 17-liter lighted reservoirs with ConSeal™ Tint Technology selectively hide surgical fluids.
- Clog-free manifold with integrated high-flow filter protects against vacuum system moisture ingress.
- Simple and interactive, high-definition touchscreen display.
- Reusable reservoirs and recyclable manifolds reduce OR red bag waste by up to 70%



# Zimmer IntelliCart System Waste Management

## Specifications

### Dimensions

**Height:** 55 in (139.7 cm)

**Width:** 23.5 in (59.6 cm)

**Depth:** 23 inches (58.4 cm)

### Display

LCD, Backlit touch screen display

### Suction/Exhaust

Closed suction system with variable suction levels based on user selection providing regulated vacuum from 0-600 mmHg using AC-powered twin piston vacuum pump passing through single-use manifolds and into canisters.

### IV Pole

Carts available with powered or manual IV poles, capable of 12,000mL maximum fluid capacity

### Safety Alert

Audible and visual display for early warning of overflow protection system and vacuum start-up alert

### Canister Exterior Graduations

Graduations marked on the canister exterior provide an approximate visual indication of contents in 100mL increments.

