



The **Zoll AED Pro** is designed with advanced capabilities and supports Basic Life Support as well as Advanced Life Support. The Zoll AED Pro features a high-resolution LCD display that allows responders to visualize a patient's ECG while utilizing a 3-lead monitoring cable. It has rubber over-molds that are built to withstand harsh environments, and utilizes the unique one-piece CPR-D padz for easy placement.

## Features

- High-quality display.
- ECG monitoring with 3-lead cable.
- Manual override options to satisfy Advanced Life Support professionals.
- Internal memory data storage and IrDA or USB data transfer.
- Feedback on CPR chest compressions utilizing the CPR-D Padz.



# Specifications

## Dimensions

**Height:** 3" (7.6 cm)  
**Width:** 9.2" (24.1 cm)  
**Depth:** 9.4" (23.9 cm)  
**Weight:** 6.5 lbs (2.9 kg)

## Defibrillator

**Waveform:** ZOLL Rectilinear Biphasic waveform.  
**Energy:** 50, 70, 85 (pediatric), 120, 150, 200 joules (adult) preconfigured.  
**Charge Time:** Less than 10 seconds with a new fully charged battery. Depleted batteries will result in a longer defibrillator charge time.  
**Energy Display:** Monitor display indicates selected energy in manual mode only.  
**AED Function:** Auto analyze and charge x3 with programmable auto energy level selection, screen prompts, and voice prompts.  
**Electrodes:** CPR-D padz, stat padz II, and Pediatric pedi padz II packaged in pairs. CPR-D padz one-piece electrode provides feedback on rate and depth of chest compressions.

## Display

**Screen Type:** High-resolution LCD display.  
**Screen Size:** 3.02 inches W x 2.27 inches H.  
**Sweep Speed:** 25 mm/sec.  
**Viewing Time:** 3 seconds.  
**Displayed Information:** Heart rate, ECG waveform, text prompts, CPR bar graph, battery gauge, elapsed time, number of shocks delivered.

## ECG Monitoring

**Patient Connection:** 3-lead ECG cable.  
**Input Selection:** Fully defibrillator protected.  
**ECG Size:** Automatic Control (AUD) gain.  
**Heart Rate:** 30-300 bpm  $\pm$ 5%.

## Battery

**Type:** Rechargeable, sealed lead acid.  
**Recharge Time:** 4 hours or less with a ZOLL PowerCharger 4x4 or 1x1 Battery Charger.  
**Operating Time:** For a new, fully charged battery pack at 20° C: 170 defibrillator discharges at maximum energy, or 6 hours continuous ECG monitoring.  
**Type:** Sealed, disposable lithium manganese.  
**Operating Time:** 300 defibrillator charges at maximum energy or 13 hours continuous ECG monitoring.

