



The **Zoll SurePower Charging Station** is a 4-bay, Zoll battery charging and testing system. The SurePower 4-bay charger allows you to charge, test, and recalibrate four Zoll rechargeable defibrillator batteries. The charging station works with both Zoll Lithium-ion and leads acid batteries. The Zoll charger offers three different charging methods: AutoTest, QuickCharge, and ManualTest. During these charging methods, the charger measures the battery's capacity to make sure it can meet the needs required for a Zoll Defibrillator to run in case of an emergency. It can charge batteries such as; SurePower, SurePower II, 4410, 4410 Smart, 4410 SmartReady, and XL defibrillator batteries. The SurePower batteries are for the Zoll R series, E series, and AED Pro defibrillators. While the SurePower II batteries are for the Zoll X Series, Propaq MD defibrillators, and Propaq M Monitor.



# Specifications

## Dimensions

Height: 6 in (15.2 cm)  
Width: 11 in (27.9 cm)  
Length: 16 in (40.6 cm)  
Weight: 15 lbs (6.8 kg) without batteries

## Power

Power requirements  
100-240 Vrms, 50/60 Hz  
120 Vrms, 400Hz  
Power Consumption: 200 VA maximum

## Temperature

Operation: 32° F to 122°F (0°C to 50°C)  
Optimal battery charging: 59°F to 95°F (15°C to 35°C)  
Storage and Transportation: -40°F to 158°F (-40°C to 70°C)  
Humidity: 5% to 95% (relative humidity - non-condensing)

