



The **Amsco Eagle 3000 Series** Sterilizers feature an advanced microcomputer control system, providing the latest standard for cycle setup, selection and monitoring. Once settings are made and the cycle started, microcomputers accurately monitor and control system operations. The Control Column includes an easy-to-read vacuum fluorescent display and an integral thermal printer.

## Features Two Configurations

### **Graviy**

For sterilizing liquids, in flasks with vented closures, at 212° to 254°F (100° to 122°C), and also heat and moisture stable goods at 212° to 280°F (100° to 138°C).

### **Prevacuum**

For efficient, high-volume processing of heat- and moisture-stable materials, such as fabrics, wrapped hard goods, glassware or animal bedding at 270° to 280°F (132° - 138°C) through prevacuum cycle. Prevacuum also provides gravity-type cycles.



# Specifications

## Single Door Product Configurations

**Amsco Eagle 3031:** 24 x 36 x 36 Gravity

**Amsco Eagle 3041:** 24 x 36 x 48 Gravity

**Amsco Eagle 3051:** 24 x 36 x 60 Gravity

**Amsco Eagle 3033:** 24 x 24 x 36 Prevacuum

**Amsco Eagle 3043:** 24 x 36 x 48 Prevacuum

**Amsco Eagle 3053:** 24 x 36 x 60 Prevacuum

## Door Hinge

On Right Side

On Left Side

## Mounting

Cabinet Enclosed / Freestanding

Recessed

