



The **Amsco Eagle 3000 Series** Sterilizers feature an advanced microcomputer control system, providing the latest standard for cycle setup, selection and monitoring. Once settings are made and the cycle started, microcomputers accurately monitor and control system operations. The Control Column includes an easy-to-read vacuum fluorescent display and an integral thermal printer.

Features Two Configurations

Graviy

For sterilizing liquids, in flasks with vented closures, at 212° to 254°F (100° to 122°C), and also heat and moisture stable goods at 212° to 280°F (100° to 138°C).

Prevacuum

For efficient, high-volume processing of heat- and moisture-stable materials, such as fabrics, wrapped hard goods, glassware or animal bedding at 270° to 280°F (132° - 138°C) through prevacuum cycle. Prevacuum also provides gravity-type cycles.



Specifications

Single Door Product Configurations

Amsco Eagle 3031: 24 x 36 x 36 Gravity

Amsco Eagle 3041: 24 x 36 x 48 Gravity

Amsco Eagle 3051: 24 x 36 x 60 Gravity

Amsco Eagle 3033: 24 x 24 x 36 Prevacuum

Amsco Eagle 3043: 24 x 36 x 48 Prevacuum

Amsco Eagle 3053: 24 x 36 x 60 Prevacuum

Door Hinge

On Right Side

On Left Side

Mounting

Cabinet Enclosed / Freestanding

Recessed

