



The Aaron 3250 New Digital System by Bovie®. The Aaron 3250 is a 300 watt multipurpose electrosurgical generator for use in the modern operating room and surgicenter. It features both monopolar and bipolar functions to satisfy all of your surgical demands with safety, flexibility, reliability, and convenience.

## Features

- Nine Presets
- Large illuminated digital displays
- Six Output Modes
- Self-test circuits
- Audible tones
- Discreet outputs
- Isolated circuitry
- Digital Error Detection



# Specifications

## Dimensions

Width: 31.1 cm (12.3")  
Height: 15.3 cm (6.0")  
Depth: 41.3 cm (16.3")  
Weight: < 8.8 kg (< 19 lbs.)

## Output Power

Cut I & II: 300 watts  
Blend: 200 watts 10 Steps  
Coagulation (Pinpoint): 120 watts  
Fulguration (Spray): 80 watts  
Bipolar: 80 watts

## Line Voltage

100-240 VAC  $\pm$  10%  
Output Frequency 490 kHz  $\pm$  5  
Line Frequency 50-60 Hz

