



The basic function of a 3-button handpiece is to allow the user to control power up/down and to activate the high-frequency desiccator. The Aaron 940 has a technological advantage in the handpiece: a superbly engineered collet that is gold plated for superior conductivity. The collet is attached to our circuit board and the entire handpiece is encapsulated and formed in a high tech polymer that combines comfort for you with a sure grip. The buttons on the handpiece have been designed for superior tactile feel and comfort. Finally, it's autoclavable so you can have sterility when you need it, and an easy-to-clean and durable handpiece all the way around.

Features

- More power with 40 watts
- More accurate bi-polar increase and decrease in 1/10th watt increments from zero to ten watts
- More accurate power with Universal Power Control — 100 - 240 VAC
- More user-friendly with features that can save you time and money
- Less to do with fewer steps to accomplish the same procedures
- Less to buy — fewer accessories needed for additional functionality
- Foot control activation is available and can help you keep the handpiece steady



Specifications

Dimensions

Height: 15.9 cm (6.3")
Width: 22.2 cm (9")
Depth: 11.4 cm (4.5")
Weight: < 1.4 kg (< 3 lbs.)

Output Power

To 40 watts (0-10 watts in
1/10th watt increments on
both mono and bipolar)

Line Voltage

Universal Power Supply
A940: 100 to 240 VAC

Output Frequency

550 kHz

Line Frequency

50 - 60 Hz