



The first High-Frequency Desiccator with Cut. One half of the unit is a high-frequency desiccator, The other half is a Bovie Generator, giving 60 watts of cutting power and coag capabilities – plenty of power for LLETZ/LEEP or most any cutting procedure performed in the doctor's office.

## Features

- Versatility - High Frequency Desiccator with 60 watts of Cutting Capabilities
- Unique - Two machines in one convenient package
- Simplicity
- Accuracy - Adjustable in 1/10" watt increments in the first 10 watts, for the most accurate output adjustable
- Cost Effective - No grounding pad required when working in the desiccation mode.

## Bovie Generator

- 60 Watts – Cut, Blend, Coagulation
- Isolated Power
- 6 Presets – For The Doctor's Favorite Settings



# Specifications

## Dimensions

Width: 26 cm (10.25")  
Height: 17.8 cm (7")  
Depth: 11.4 cm (4.5")  
Weight: < 3.7 kg (< 8 lbs.)

## Generator

Cut: 60 watts  
Blend: 60 watts  
Coagulation: 60 watts

## High Frequency Desiccator

Fulguration: 35 watts  
Bipolar: 30 watts

## Line Voltage

A950™: 120 VAC ± 10%  
A950™-220: 230-240 VAC

## Output Frequency

350-800 kHz

## Line Frequency

50 - 60 Hz

