



The Alaris® 8000 Point-of-Care (PC) Unit is the core of the Alaris® System and provides a common user interface for programming infusion and monitoring modules. This helps to reduce complexity at the point of care. The Medley Medication Safety System is a modular infusion and monitoring system designed to provide SpO2 monitoring capabilities and accurate, automated infusion of a broad range of intravascular fluids, medications, and blood products.

Customize this unit by adding any of the modules below

- Alaris 8100 Pump Infusion Pump
- Alaris 8110 Syringe Module Infusion Pump
- Alaris 8120 PCA Module Infusion Pump
- Alaris 8200 SpO2 Module

