



CMP Warming Cabinets maintain ideal temperatures for blankets, solution bottles and bags. CMP has been building high quality warming cabinets for over 50 years and has become the industry standard for hospital and emergency room blanket and solution storage. CMP stainless steel Warming Cabinets can be found in almost every hospital in the U.S. and can be configured for your unique requirements

Features

- Single or Double Compartment in Various Sizes
- Stainless Steel Construction
- Microprocessor Control System
- Digital Display Indicators
- Fahrenheit or Celsius Temperature Scale
- U.L. Listed



Specifications

Dimensions

Height: 37.5 in (95.3 cm)
Width: 17 in (43.2 cm)
Depth: 17 in (43.2 cm)
Weight: empty: 131 lbs. (59.4 kg) full: 147 lbs. (67kg)

Fluid Temperature Range – Manual

4°C – 42°C (39.2°F – 107.6°F)

Patient Temperature Set Point Range

30°C – 40°C (86°F – 104°F)

Reservoir capacity

2 gallons (7.5 liters)

Flow Rate

42 gph (115v); 36 gph (230v)

Heater

800 Watts

Probe Type

400 series

