



The Malis CMC-III System offers microprocessor control of an exclusive generated waveform.

Features

- High Power and Micro Bipolar Cutting. Advantages of high power as well as micro bipolar cutting of a variety of tissues, with cutting similar to mono polar
- Unique Patented Waveform. Provides efficient, gentle coagulation of tissues whether bathed in body fluids or irrigation
- Dual Foot Pedal Capabilities. MALIS CMC-III Bipolar System provides an electric foot pedal with cut and coagulation capabilities
- Precise Fluid Delivery. The easy to use MALIS Irrigation Module 1000 Bipolar Electrosurgery Irrigator uses a peristaltic pump to deliver a consistent and repeatable rate of fluid to meet a wide range of surgeons needs

