



The Codman Synergy Malis features a unique patented waveform with reduced charring/sticking and extended forcep life.

Features

- Unique patented waveform.
- Reduced charring/sticking.
- Extended forcep life.
- Ultra low output impedance.
- Smooth progressive coagulation.
- 0.045-50 Watts.
- Reliable durable unit.
- Wide range of forceps and accessories.
- Tabletop or cart mount design.
- Portable unit.



Specifications

Dimensions

Height: 13.5 cm

Width: 26 cm

Depth: 37.5 cm

Weight: 11 Lbs.

General

OutPut Frequency: 1000 kHz

Output Power: 50 W

