



With a multitude of applications in such diverse specialties as dermatology, plastic surgery, gynecology, dental and oral surgery, ophthalmology, otolaryngology, proctology, and urology, the Conmed Hyfrecator Plus 7-797 lets you perform a wide variety of desiccation, fulguration, and coagulation office procedures.

The ConMed Hyfrecator Plus 7-797 includes a foot control pencil with 10 foot (3.0 m) cord, foot pedal, and sterile sheaths. The Hyfrecator Plus 7-797's electronic circuitry provides a near linear power output. It provides precision whether it's being used at full power for broad surface coagulation or low output for delicate facial procedures.

Features

- Monopolar and bipolar outlets
- Banana jack outlet
- Small compact 10 x 5 1/2 x 4
- Hand pencil



Specifications

Dimensions

Depth: 4 in (102 mm)
Width: 8.75 in (222 mm)
Height: 7.5 in (190 mm)
Weight: 6 lbs (2.7 kg)

Environmental Conditions for Transport & Storage

Temperature Range: -20 to + 60 Degrees C
Humidity: 0% to 97%

Output Power

Maximum Power: 36 Watts (500 Ohm Load)
Bipolar: 30 +/- 6 Watts (500 Ohm Load)
High: 30 +/- 6 Watts (1000 Ohm Load)
Low: 15 + 2/6 Watts (500 Ohm Load)

Open Circuit Voltage

bipolar: 3000 Vp-p Max.
High: 8000 Vp, 10,000 Vp-p Max.
Low: 4000 Vp-p Max

Leakage

60 hz: 50 microamps max, with third wire of line open.

Frequency

500 Khz at 31.25 KHz rate