



The **Datex Ohmeda Modulus SE Anesthesia Machine** Provides a set of medical gasses to a breathing system. The Modulus SE can be equipped with either the 7800 ventilator or the 7900 Ventilator. The System offers a large selection of gases, vaporizers, ventilators and monitors that allows you to configure the system for what is needed for each situation. The Modulus SE anesthesia Machine Comes with 3 standard drawers and an overhead shelf for storage.

## Features

- Ascending bellows ventilator (7800 or SmartVent 7900). Modes CMV, PCV (with SmartVent), Manual/Spontaneous.
- Mechanical/pneumatic flowmeters (needle valves and glass flowtubes) and vaporizers.
- Traditional size, configuration of absorber head and higher-volume breathing circuit. Loose absorbent granules or canisters.
- Removable vaporizers. Can be used with non-rebreathing circuits. Opened or closed scavenger interface.
- Monitoring Included: Volume Pressure and Oxygen Inspired.



# Specifications

## **Dimensions**

Height: 170 cm (66.9")  
Depth: 79.4 cm (31")  
Width: 79.4 cm (31")  
weight: 139.3 kg (309.5 lbs)

## **Drawers and shelf**

Drawer: 47.6 x 34.3 x 14 cm (18.8" x 13.5" x 5.5")  
Fixed Shelf: 47.6 x 36.8 cm (18.8" x 14.5")  
Slide Shelf: 46.5 x 36.8 cm (18.3" x 14.5")

## **Absorber and Post**

Absorber Arm Length: 30.5 cm (12")  
Vertical Adj (above floor): 39.4 to 71.1 cm (15.5" to 28")