

# Datex Ohmeda Tec 5 Vaporizer



## Sevoflurane



### Features

- It offers significant performance and convenience benefits, coupled with important safety features, to improve the clinical management of volatile
- Anesthetic agents in continuous flow techniques.
- Options
  - Forane
  - Ethrane
  - Halothane
  - Sevoflurane



### Specifications

#### Dimensions

Height:	237mm
Width:	114 mm
Depth:	97 mm
Weight:	7.0 k

#### Calibration:

The dials of all Tec 5 Vaporizers are calibrated to 5% v/v in steps of 0.2% up to 1%, and in steps of 0.5% from 1% to 5%

#### Operating Temperature:

18-35C

#### Resistance to gas flow:

“OFF” Isolated- no resistance from vaporizer  
“ON” 5 L/min O<sub>2</sub> 15-20 cmH<sub>2</sub>O

#### Liquid capacity:

To charge dry vaporizers: 300mL (nominal)  
Retained by wick system: 75mL (nominal)