

Datex Ohmeda Tec Plus Vaporizer

Desflurane



Features

- Designed specifically for the delivery of Desflurane anesthetic agent
- Heats Desflurane to maintain temperature and vapor pressure for consistent output

Specifications

Dimensions

Height:	9.9"
Width:	9.3"
Depth:	4.3"
Weight:	20 lbs

Concentration range:	1 to 18%
Calibrated flow range:	0.2 to 10 L/min (at 21 degrees Celsius)
Max. permissible flow:	20 L/min
Liquid capacity:	375 mL working volume 250 mL refill indicator 50 mL reserve
Voltage requirement:	either 220 to 240 V or 100 to 120 V a.c. 50/60 Hz
Max. permissible flow:	Duracel MN/1604

