



The ERBE EIP 2 provides a pulsatile jet of sterile water for rinsing substances such as blood and fecal matter from the procedural site and instruments.

Features

- Works in open endoscopic procedures,
- Uses only disposable ERBEFLO accessories for optimum time savings,
- Pre-set buttons offer customized flow rates,
- Clear flow rate display,
- High volume irrigation,



Specifications

Dimensions (not including pump head)

Width: 8.1" (20.5 cm)
Height: 4.9" (12.5 cm)
Depth: 6.5" (170 cm)
Weight: 6.6 lbs. (3 kg)

Power Connection

Rated supply voltage 120 V \pm 10%
Line current 400 mA
Power fuses T 0.4 A
Rated supply frequency 50/60 Hz
Power input in standby mode 3 Watt
Power input with max. HF output 20 Watt
Terminal for grounding (potential equalization) Yes

Operating

Flow rate range 20% (100 ml/min) to 100% (500 ml/min) \pm 50 ml/min
Pre-sets and adjustment increments 3 pre-set flow rates and 1% increment adjustments
Continuous Yes, but for only a maximum time of 20 seconds