



Designed especially for the GI Endoscopist. The units are particularly useful during Colonoscopy with active hemorrhage since it allows the delivery of a pressurized lavage to improve visualization during difficult procedures.

Features

- Front of console has peristaltic pump with plastic tube guide and enclosure
- Variable flow dial
- Pneumatic footswitch connector; footswitch not included
- Power on/off rocker switch with green indicator light
- Universal power cord



Specifications

Dimensions

Width: 5.5 in
Depth: 10
Height: 4.5

Electrical

Volts 95-135 V, 185-250V
Hertz 50-60 Hz

Pump

Flow Rate: 150-650 ml/min
Pressure: 2-40 psi
Reservoir: 500 ml Removable Bottle

