



The APC 300 is the 1st generation Gastrointestinal Intervention system. The Argon Plasma System redefines interventional treatment options in flexible endoscopic intervention, both gastrointestinal and tracheobronchial. Cut control - for automatic perfect cuts

Features

Cutting

- High cut with specialized arc regulation for underwater cuts or in fatty structures
- Endo cut - the cutting mode for endoscopy
- Mono- und bipolar cutting (Bipolar Cut)
- Power Peak System (PPS) for intelligent support of initial incision and cutting

Coagulation

- Precisely controlled coagulation and no carbonization with Soft Coag
- Forced Coag - for situations where speed is essential
- Spray Coag - for efficient surface coagulation
- Argon-Plasma Coagulation unit ERBE APC 300 hookup as standard equipment
- Mono- und bipolar coagulation (Bipolar Coag)
- Twin Coag - for comfortable coagulation on two channels:connector and individual settings for two Coag instruments
- Auto Start and Auto Stop - for coagulation support

Handling

- Maximum user-friendly operation via 11 individually configurable programs

Safety

- Nessesy the one-of-a-kind Neutral Electrode Safety System
- Double-safe cutting with Cut Control and Bipolar Cut



Specifications

Dimensions

Width: 16,14 in (410mm)

Height: 4.1 in (105 mm)

Depth: 15 in (380 mm)

Weight: 15 lbs. (6.8kg)

Power

Power input 230 V

Frequency 50-60 Hz

Power 28 W

Ambient Conditions for transport and Storage

Temperature -40 °C to +70 °C

relative Humidity: 10% to 95%

Ambient Conditions for Operation of Unit

Temperature: +10°C to +40°C

Relative Humidity: 5% to 80%, non-condensing

Cutting

Cut control - for automatic perfect cuts

High cut with specialized arc regulation for underwater cuts or in fatty structures

Endo cut - the cutting mode for endoscopy

Mono- und bipolar cutting (Bipolar Cut)

Power Peak System (PPS) for intelligent support of initial incision and cutting

Coagulation

Handling

Maximum user-friendly operation via 11 individually configurable programs

Safety

Nessy the one-of-a-kind Neutral Electrode Safety System

Double-safe cutting with Cut Control and Bipolar Cut