

Datex Ohmeda Tec 7 Vaporizer

Sevoflurane



Features:

- Delivers Isoflurane, Sevoflurane, Enflurane, and Halothane anesthetic agents
- Provides output consistent with the dial setting throughout the clinical flow range
- Features an easy-to-turn dial to help control agent delivery
- Offers two filling options: Easy-fil* or Quik-Fil* (Sevoflurane only)
- Simplifies agent filling and helps minimize agent leaks while filling

Specifications:

- Tec 7 Vaporizers are designed to provide consistent output throughout the clinical flow range from 200 mL/min to 15 L/min.
- Models are available for Isoflurane, Sevoflurane, Enflurane, and Halothane.
- Tec 7 Vaporizers are equipped with a large diameter control dial, which incorporates fine graduations of 0.2% between 0 and 1%, and 0.5% from 1% to 8%. The dial for Sevoflurane is marked in steps of 0.2% up to 1% v/v, and in steps of 1% between 1% and 8%.
- The easy-turning dials and small graduations on the Tec 7 Vaporizer help you fine tune anesthetic delivery over the range of dial settings and flow rates.
- To help minimize filling frequency, the Tec 7 Vaporizer accommodates 225 mL of anesthetic agent.
- The Tec 7 Vaporizer can be mounted on GE Healthcare anesthesia delivery systems equipped with the Selectatec* manifold without tools and without taking the anesthesia system out of service.
- An innovative non-spill system in the Tec 7 Vaporizer limits movement of liquid agent if the vaporizer is tilted or inverted, helping to protect internal components.