



The **HOTLINE®** blood and fluid warmer offers a safe, reliable, and proven technology to consistently warm and deliver blood and fluids to patients undergoing surgical procedures. Fluid warming is a method of raising the temperature of fluids administered to a patient to maintain normal body temperature and prevent hypothermia.

## Features

- 42°C recirculating solution temperature set point ensures the patient receives blood or fluids at a temperature able to maintain normothermia
- The **HOTLINE®** triple lumen disposable design eliminates patient line cool-down, resulting in warmer fluids reaching the patient
- Keeps blood and fluids warm between 37°-42°C
- Easy to use and simple, one-step disposable setup Outperforms hotplate technology at routine flow rates

