



The Aspect 100-4T provides a very large radiolucent imaging area for easy C-arm access and greater anatomical coverage. A high-quality, 24" wide, carbon fiber tabletop and 2" thick, specially constructed tabletop pad provide for minimal x-ray attenuation and maximum patient comfort.

Features

- Cantilevered, carbon fiber tabletop with low x-ray attenuation
- Choice of tabletop dimensions:
 - 23.5/30" x 80" (60/76 x 203 cm),
 - 18/30" x 80" (46/76 x 203 cm)
- Unobstructed imaging area: 72" (182.9 cm)
- Motorized Elevation 27.9" – 37.7" (71 – 95.8 cm)
- 4-way manual floating tabletop with panning hand control
- Longitudinal tabletop travel: 40" (101.6 cm)
- Transverse tabletop travel: 10" (25.4 cm)
- Locking swivel casters
- Patient weight capacity: 450 lbs. (204 kg)
- Table weight: 490 lbs. (222 kg)
- 120 VAC or 230 VAC operation (must specify on order)
- Battery back-up operation



Specifications

| | |
|------------------------------|---|
| General | Carbon fiber tabletop with low x-ray attenuation 2" thick tabletop pad with memory foam construction Rail-mounted panning hand control for elevation and 4-way float functions 5" (12.57 cm) locking swivel casters |
| Table Dimensions | Choice of two tabletops: Contoured, 23.5" / 30" x 80" (59 cm / 76 cm x 203 cm) Contoured, 18" / 30" x 80" (46 cm / 76 cm x 203 cm) |
| Table Movements | 72" (182.9 cm) unobstructed imaging area at maximum longitudinal float Motorized Elevation 27.9" – 37.7" (71 – 95.8 cm) 4-way manual floating table to with panning hand control Longitudinal tablet float 40" (101.6 cm) Transverse tabletop float 10" (25.4 cm) |
| Electrical | 110 VAC, 60 Hz operation (230 VAC operation optional). ≤10 amp @ 110V. Battery back-up operation |
| Table Weight Capacity | 450 lbs (204 kg) |
| Table Weight | 490 lbs (222 kg) |
| Accessories included | Tabletop Pad & Patient Restraint Straps |
| Optional Accessories | Table mounted radiation shield Radiolucent arm boards Catheter Tray Extension |

