



Multi-Purpose Mobile Urological Imaging Table with motorized Elevation, Trendelenburg Tilt Longitudinal and Transverse Tabletop Motions and carbon fiber tabletop with low x-ray attenuation.

Features

- Designed for endourology procedures
- 4 movement
- 10" longitudinal slide
- 250-525 lb capacity
- 33.3"-43.1" high-low
- trendelenburg +/-15°
- transverse slide +/- 3.5"



Specifications

Table dimensions and components

- Carbon fiber tabletop with low x-ray attenuation- 2" (5 cm) thick tabletop pad and tabletop extension pad with memoryfoam construction
- Tabletop: 28"W x 46"L (71 x 117 cm)
- Radiolucent, lightweight tabletop extension: 24"W x 34"L (61 x 86 cm)
- Built-in footswitch storage tray
- Removable radiographic cassette tray, loads from either side of table
- Standard US side-rails for accessory attachments on main tabletop
- Metal-free perineal end of tabletop for unobstructed fluoroscopy access to edge of table

Table movements

- Locking swivel casters, 5" (12.7 cm), and built-in transport handle
- Motorized motions via hand control pendant and foot switch control
- Elevation: 31.3" – 41.1" (79.4 – 104.4 cm)
- Trendelenburg Tilt: ± 15°
- Longitudinal tabletop travel: 10" (25.4 cm)
- Transverse tabletop travel: ± 3.5", 7" total (± 8.9 cm, 17.8 cm total)

Imaging area

- Imaging area, main tabletop: 25"W x 20"L (63.5 x 50.8 cm)
- Imaging area, tabletop extension: 24"W x 34"L (61 x 86 cm)
- Total radiolucent imaging area length with extension: 54" (137 cm)

Patient weight capacity

- With patient upper body on main tabletop: 525 lbs.(238 kg)
- With patient upper body on extension: 400 lbs. (181 kg)
- Load capacity on extension only: 200 lbs. (90.7 kg)

Electrical

- 110 VAC, 60 Hz operation (230 VAC operation optional). ≤10 amp @ 110V

Table & Extension Weight

- Table: 670 lbs. (304 kg)
- Tabletop Extension: 16 lbs. (7.25 kg)

Accessories included

- Tabletop Pad, 2" thick
- Tabletop Extension, with 2" thick pad
 - Hand control pendant & footswitch control
 - Removable Drain Bag Support Hoop
 - Disposable drain bag and hose ass'y, box of 20
 - Removable Radiation Shield Flap for perineal end of table
 - Patient Restraint Straps (2)
 - Armboard, rail mounted, with pad and strap
 - Disposable covers for footswitch, box of 50

