

Specifications

Mechanical and Electrical Specifications

Carbon fiber tabletop with low x-ray attenuation

Tabletop = 24" x 46" (61 x 117 cm)

22" W x 32"L (56 x 81 cm) unobstructed imaging area

Motorized tabletop elevation with foot control:

- Elevation 27" – 39"
- (68.5 – 99cm) with tabletop pad
- Locking swivel casters
- Table weight ≈ 375 lbs. (≈170 kg)
- 120 VAC or 230 VAC operation 10 amp max. required at 110V (must specify on order)

Patient Weight Capacity

450 lb. (204 kg) patient weight capacity

