



The Aspect ISR G2 table for vascular and interventional procedures is IDI's most advanced C-arm imaging table. It delivers extended positioning and imaging capabilities at an affordable price. The ISR has 32" of tabletop float for very large anatomical imaging coverage without C-arm movement.

Features

- 4-way Manual-Float Tabletop
- Motorized Elevation
- Iso-Roll
- Trendelenburg tilt
- Reverse Trendelenburg Tilt
- Choice of tabletops



Specifications

Table Movements

Elevate and Tilt functions are fully motor driven - no hydraulics!

Elevation: 32" - 40" (81-102cm), including 2" (5cm) tabletop pad.

Trendelenburg Tilt: $\pm 12^\circ$ with automatic "stop-at-level" positioning.

Lateral roll: $\pm 12^\circ$, isocentric.

4-way manual floating tabletop with panning hand control.

Longitudinal tabletop travel: 32" (81.3 cm).

Transverse tabletop travel: 8" (20.3 cm).

Automatic locking of longitudinal float when table is not level (head up or head down). Transverse float is available regardless of tilt or roll position.

Dual-wheel 6 caster system with 3-position locking system: total unlock, total lock, and steering lock (for easy transport).

Imaging Area

Imaging area at full longitudinal extension = 73" (185cm) with optional headrest attached and 61" (155 cm) without optional headrest attached.

Table Dimensions

Standard tabletop dimensions: 19.5"/22"/24" X 84"

(49.5cm/55.8cm/60.96cm x 213.3 cm). #X100-2492

Optional tabletop dimensions: 22"/24" X 84"

(55.8cm/60.9cm x 213.3 cm). #X100-2426

Weight

Patient weight capacity: 450 lbs. (204 kg) with full table functionality at max. patient weight.

Table weight: 605 lbs. (274 kg).

Electrical

Full operation on either AC or internal battery power.