



Designed for hospital induction, shock treatment, and recovery rooms, the LUCEA 40 offers an illumination of 40,000 lux. Its 26 LEDs achieve far higher shadow dilution than other lights in its category.

Features

- As well as mobile lights, LUCEA 40 lights are available in ceiling and wall-mounted versions which free the floor space and facilitate cleaning
- Having a design with no openings also contributes to improved hygiene
- The LUCEA 40 is also designed to comply with the most stringent standards (IEC, UL, etc)



Specifications

General

Size of lighthead (mm / inches) 337 x 214 / 13.3 x 8.4

Size of lighthead (kg / ilbs) 19.1 / 4.1

Illumination at 1m / 39.4 inches (lux) 40,000

Field diameter 1m / 39.4 inches (cm/inch) 24 / 9.4

Color temperature (k) 4,500

Color rendering index (CRI) (RA) 95

Rated power 55

Lighthead power consumption (watt) 20

Nominal voltage 45