



The PowerLED 300/300 surgical lights provide a combination of two PowerLED 300 surgical lights. The PowerLED 300 surgical lights are 17 cm light patch for deep and narrow cavities, heat-free and 110,000 Lux

**2nd Generation of LEDs**

- Superior luminous efficiency
- Optimized service life up to 40,000 hours
- Low energy consumption

**Flux Stability Program:**

- Patented LED luminous flux management
- Constant illumination even during extended operations

**A new ergonomic to work with:**

- Slim profile for easy handling
- A compact and featherweight light head for perfect positioning
- Central handle available as an option



## Specifications

**Optical and photometric  
characteristics: PowerLED  
300 (ea)**

Illumination (lux) 110,000  
Field diameter\* (cm/inch) 17 / 6.7  
Color temperature (k) 3,800  
Color rendering index (CRI) 95  
Radiant energy (average) (mw/m<sup>2</sup>.lx) 3.6  
LED service life (hours) 40,000

**Ambient light  
characteristics: PowerLED  
300 (ea)**

Illumination (lux) 500

