

The PowerLED 300 surgical lights are 17 cm light patch for deep and narrow cavities, heat-free and 110,000 Lux
2nd generation of LEDs:

- Superior luminous efficiency
- Optimized service life up to 40,000 hours
- Low energy consumption

- Patented LED luminous flux management
- Constant illumination even during extended operations

- Slim profile for easy handling
- A compact and featherweight light head for perfect positioning
- Central handle available as an option



Specifications

General

Illumination (lux): 110,000
Field diameter* (cm/inch): 17 / 6.7
Color temperature (k): 3,800
Color rendering index (CRI): 95
Radiant energy (average) (mw/m2.lx): 3.6
LED service life (hours): 40,000
Illumination (lux): 500