



The **Maquet MODUTEC** ceiling supply units are the most compact configurations of strong rigid arms, hygienic equipment managers and fully functional columns developed to match any clinical application need within care areas and operating theaters. The elegant MODUTEC design offers many incremental options to optimize both clinical staff and patient environments.

Features

- MODUTEC full-color LEDs ambient light embedded in the distributor offers the ultimate user experience.
- This “LED inside” option enhances the multimedia viewing experience by improving perceived picture details on OR screens.
- Conceptualized to remove dazzling light effect, it is effectively used during mini-invasive interventions for staff precise viewing of equipment on shelves.
- In care areas, the environment becomes safer and more comfortable.



Specifications

MODUTEC DISTRIBUTION MODULES CAPACITY FOR OUTLETS AND CONDUITS

Module Ranges	Max. Conduits	Distribution Modules	Length (inch/mm)	Number of Panels	Electrical Outlets Per Panel	Gas Outlets Per Panel	Cable Channel Outputs Per Module
TONIC	8	1001	39.4/1000	1	6	7	7 opt.
		1201	47.2/1200	1	9	11	11 opt.
		1401	55.1/1400	1	9	11	11 opt.
		1601	63/1600	1	12	15	15 opt.
		1801	70.9/1800	1	12	15	15 opt.
		2001	78.7/2000	1	12	15	15 opt.
		2201	86.6/2200	1	12	15	15 opt.
QUADRA	8	402	15.7/400	2	2	3	4 incl.
		802	31.5/800	2	6	7	6 incl.
		1202	47.2/1200	2	9	11	10 incl.
		1602	70/1600	2	12	15	12 incl.
MOBILI	12	403	15.7/400	3	2	3	3 opt.
		803	31.5/800	3	6	7	7 opt.
		1203	47.2/1200	3	9	11	11 opt.
NEST	16	405	15.7/400	5	2	3	4 opt.
		805	31.5/800	5	6	7	6 opt.
		1205	47.2/1200	5	9	11	10 opt.
RACK	16	500	27.6/700	7	1, 2 or 3	1, 2	1, 2 opt.
		800	39.4/1000	7	1, 2 or 6	1, 2	1, 2 opt.

MODUTEC EXTENSION MODULES CAPACITY FOR OUTLETS AND CONDUITS

Extension Ranges	Max. Conduits	Distribution Modules	Length (inch/mm)	Number of Panels	Electrical Outlets Per Panel	Gas Outlets Per Panel	Cable Channel Outputs Per Module
TONIC	8	401	15.7/400	1	2	3	3 opt.
		801	31.5/800	1	6	7	7 opt.
		1201	47.2/1200	1	9	11	11 opt.
QUADRA	8	402	15.7/400	2	2	3	4 incl.
		802	31.5/800	2	6	7	6 incl.
		1202	47.2/1200	2	9	11	10 incl.

GROSS LOADINGS ON MODUTEC EXTENSION ARMS

Range of Arm	6.0	8.0	9.0	10.0	12.0	6.6	6.8	9.6	6.10	9.8	9.9	9.10	12.6	12.8	12.9	12.12
Reach (in)	23.6	31.5	35.4	39.3	47.2	47.2	55.1	59.1	63	67	70.9	74.8	70.9	78.7	82.7	94.5
Reach (mm)	600	800	900	1000	1200	1200	1400	1500	1600	1700	1800	1900	1800	2000	2100	2400
Light Arm Type LL (lbs)	685		509				387		308				253			
Medium Arm Type MM (lbs)	771		679		533	518		412			339		335		277	233
Strong Arm Type HM (lbs)	771		771		641						425				352	300
Heavy Arm Type HH (lbs)						771		771			760		749		638	550
PL G 12-0 Spring Arm (lbs)					66											
Power Arm Type PWR (lbs)		441		331			441		331	441		331		441		