



Featuring innovative Microstream™ CO2 technology and reliable Nellcor pulse oximeter, the units offer many of the same advantages as the capnography and pulse. Each rugged device weighs less than 2 lbs and measures approximately 8"x 3.5"x 2".

Capnography

Display Range: 0-99 mmHg; 0-13.0 Vol %; 0-13.2 kPa
CO2 sampling flow rate: 50 ml/ min
Warm-up time: 30 seconds typical
Accuracy: 0-38 mmHg ± 2 mmHg; 39-99 mmHg ± 5% of reading
Respiration rate: 0-150 breaths/min
Modes: Adult, neonatal

Pulse Oximeter

Display Range: SpO2 Saturation 0%-100%
Pulse Rate: 20-250 beats/minutes
Accuracy: Adults 70%-100% ± 2 digits, 0-69% unspecified
Neonate: 70%-100% ± 2 digits, 0-69% unspecified

