



The Midmark 405 Power Exam Table is a versatile tool in any work environment because it's so useful to a wide range of practice groups, like obstetricians, gynecologists, dermatologists, and plastic surgeons. The reason is because of its hand-controlled Power Up/Down and Back functions, adjustable stirrups that lock in 1 of 3 lateral positions, four storage drawers, foot extension shelf, paper roll holder, removable upholstery, treatment pan, four drawers and pelvic tilt.

Features

- Auto Return (HOME) Function on Hand Control
- Disable Switch, which allows operator to disable table during a procedure
- Pelvic Lift Bar allowing seat section to be positioned in Lithotomy position
- Four storage drawers on steel ball bearing glides; two drawers on foot end of table and two drawers on either right or left side of table
- Styled drawer fronts which have molded hand pulls
- Duplex receptacle, located on left side of table
- Paper roll holder which is hidden in head end of upholstered top
- Optional Foot Control, which allows operator to position table or initiate the AUTO RETURN function, without using hands



Specifications

Dimensions

Styled Upholstered Top: 27" W x 55" L
Length with foot extension: 73"
Table Width: 27"
Minimum Height: 25"
Maximum Height: 37"
Base: 18" x 35-1/4"

Exam Table Specifications

Weight: 454 lbs
Back Section Adjustment: 0° to 70°
Table Top Height Adjustment: 25" – 37"
Pelvic Lift Adjustment: Seat section adjustable to +5°.
Stirrup Extension (Maximum): 14"
Stirrup Lateral Positions: Four lateral positions
Footrest Extension (Maximum): 17-3/4"
Footrest Extension: Adjustable to +7°

Electrical Requirements

Power Requirements :115 ±10% VAC, 50/60 HZ, 12.0 amps
Power Cord: Extends 2.5 m (98.0 in.) 12/3 AWG, Type SJT 300 V cord with hospital grade grounding type plug