



Positioning patients can be a difficult task when dealing with less mobile, elderly and special needs patients. The Ritter 230, with its low height and multi positioning capabilities is designed to improve office efficiency by allowing you and your staff to easily position patients with the touch of a button—reducing the likelihood of fatigue, back and neck strain, and workplace injury.

The Ritter 230 is designed to make patients feel more comfortable and secure when they can easily get onto a table and experience smooth and steady movement during positioning. This allows patients' anxiety to be reduced and helps them feel satisfied with the level of care you provide.

Consider the Ritter 230 Power Procedures Table a long-term investment for your practice, your staff and your patients.

Features

- A power procedures table with all the basic features you need
- Powered multi-positioning capabilities, including knee-chest
- Ergonomic, intuitive standard hand control in an easy-to-find location



Specifications

Dimensions

Minimum 22.5" (57.2 cm)
Maximum 40" (101.6 cm)
Weight: 450 lb (204 kg)

Standard Features

Stirrups
Paper roll holder holds 18" x 3.5" paper roll (45.7 x 8.9 cm)

Options

Rotation: models equipped with rotation option add 3 1/2 inches to minimum and maximum height specifications
Electrical outlets: duplex, hospital grade (on each side)
Flat upholstery set, 24" wide

Electrical

115 VAC, 10 Amps, 50/60 Hz (for electrical outlet models)
115/230 VAC, 7/3.5 Amps, 50/60 Hz (for non-electrical outlet models)
115/230 VAC, 7.5/3.8 Amps, 50/60 Hz (for rotation models)

