



The Midmark 491 power table positions patients exactly how you want them at the touch of a button. The low height ability of the chair, at only 22", is great for patient getting on and off the chair, as well as positioning them where you need them during a procedure. The power controls are on both sides of the chair, making adjustment easy and convenient. An optional foot control is also available. The back is operated at the touch of a button, as well as the height adjustment and the auto return feature. The base also rotates 330° for even more positioning flexibility. The chair is extremely stable, helping both the physician and patient feel confident and safe.

Features

- Low height ability
- Optional foot control is also available
- Back, height adjustment and auto return feature are operated at the touch of a button
- Base rotates 330°
- Movements and Features: Table Up/Down, Back Up/Down, Auto Return, and Arm Rests



Specifications

General

Minimum Height: 22" (55.9cm)

Maximum Height: 31" (78.7cm)

Width (Min): 24" (61cm)

Width (Max): 26.5" (67.3cm)

Rotation: 330o

Electrical: Domestic: 115 VAC, 60Hz, 12AMP

