



This chair boasts a very low minimum height of 19". This low entry and exit height make it very easy for the patient to get on and off the chair safely and confidently. This power chair comes with a footswitch or hand control for operating movements. The Midmark 630 also comes standard with an accessory side rail, paper holder, and articulating headrest. The Midmark Ritter comes with a pull-out drain pan at the end of the seat section and pulls out heel stirrups. Not only can the Midmark 630 be used in Family and General Practice, but also in OB/GYN, Internal Medicine, Dermatology, Plastic Surgery, Colon and Rectal, General Surgery, and other specialties

Features

- Table Up/Down
- Tilt
- Foot Up/Down
- Back Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper holder
- Table paper straps
- Foot Switch or Hand Control
- Accessory Side Rails



Specifications

General

Minimum Height: 19" (48.26cm)

Maximum Height: 40" (101.5cm)

Table Length (foot extended): 72" (182.9cm)

Table Length (foot and headrest extended): 83" (210.8cm)

Width (seat): 28" (68.6cm)

Width (back): 24" - 26" (58.4-66.0 cm)

Patient Weight Capacity: 450 lbs.(204.1 kg)

Paper Roll Holder: Holds 18" x 3.5" paper roll (45.7cm - 8.9cm)

Foot or Hand Control: To power height, tilt, back, foot/leg, Stop, auto return

Electrical - Domestic: 115 VAC, 60Hz, 10AMP, UL Listed, CSA and CE Certified