



The **Midmark Ritter 75e Evolution** power procedure chair/table is an excellent choice for multiple types of procedures. the height of the chair can be adjusted from 26 inches to 42 inches. the 75E Procedure chair can achieve a trendelenburg angle of 45 degrees. This power chair comes with a foot switch for operating movements, accessory side rail, paper holder, and articulating headrest. The Midmark Ritter 75E comes with a pullout drain pan at the end of the seat section and pull out heel stirrups. The foot/leg rest section of the 75 Evolution is lockable in 4 proctology chair positions.

Features

- Movements: Table Up/Down, Tilt, Foot Up/Down, and Back Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper Holder
- Table Paper Straps
- FootSwitch
- Accessory Side Rails.



Specifications

Dimensions

Minimum Height: 26" (66 cm)
Maximum Height: 42" (106.7 cm)
width (top): 27" (68.6 cm)
Weight: 475 lbs.

Tilt

Trendelenburg to 45 Degrees

Paper Roll Holder

21" (53.3cm)

Foot Control

To power height, tilt, back & foot

Electrical

Domestic: 120 VAC, 60Hz, 12AMP, UL Listed, CSA Certified

