



This chair boasts a very low minimum height of 22.5". This makes it very easy for the patient to get on and off the chair safely and confidently. This power chair comes with an 8-way footswitch or hand control for operating movements, accessory side rail, paper holder, and articulating headrest. The Midmark Ritter 75L comes with a pull out drain pan at the end of the seat section and pull out heel stirrups. The foot/leg rest section of the 75L is lockable in 4 proctology chair positions

Features

- Movements: Table Up/Down, Tilt, Foot Up/Down, and Back Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper Holder
- Table Paper Straps
- Foot Switch or Hand Control
- Accessory Side Rails



Specifications

Dimensions

Upholstered Top: 27" (W) x 69.25" (L) (68.6 cm x 175.9 cm).

With Headrest Extended: 82.75" (210.2 cm).

Minimum height: 22.5" (57.15cm)

Maximum height: 40" (101.6 cm)

General

Patient Weight (Maximum): 325 lb (147.4 kg).

Weight of Table: 515 lb (233.6 kg).

Back Section Range: 0 (horizontal) to $85^{\circ} \pm 5^{\circ}$.

Table Top Height Range: $22.5" \pm 0.5"$ to $40" \pm 1"$ (55.9 cm \pm 1.3 cm to 101.6 \pm 2.5 cm).

Table Top Tilt Range: 0 (horizontal) to $30^{\circ} \pm 5^{\circ}$ (foot up).

Foot and Head Controls

To power height, tilt, back, foot/leg, Stop, auto return

Electrical

115 VAC nominal (126 VAC maximum - 110 VAC minimum), 60 Hz, 12 amps maximum.

Extends 70" (minimum) from table.

12 AWG / 3 conductor.